

AFZAN LAB MANUAL

Experiment No. 01

Page No.

Date 31-07-26.

1. Preparation of a True Solution of Sugar in Water.

- Requirements:
 - 1) Beaker
 - 2) Glass Rod
 - 3) Water
 - 4) Sugar
 - 5) Spoon

- Procedure:
 - 1) Take about 100ml of Water in a beaker.
 - 2) Add one spoon of sugar to the Water.
 - 3) Stir the mixture with glass rod until the sugar dissolves completely.
 - 4) Observe the solution.

- Observation:
 - 1) The sugar dissolves completely in water.
 - 2) The solution is clear and transparent.
 - 3) No particles are ~~not~~ visible.

- Result: A true solution of sugar in water was prepared successfully. The sugar dissolved completely in water, forming a clear and transparent solution. No undissolved particles were ~~not~~ visible, showing that sugar is completely soluble in water.

AFZAN LAB MANUAL

Experiment No. 02

Page No.

Date 31-07-26.

2. Preparation of Suspension of Chalk Powder (or Mud) in Water:-

- Requirements:-
 - 1) Beaker
 - 2) Water
 - 3) Chalk powder (or Mud)
 - 4) Glass Rod
 - 5) Spoon

- Procedure:-
 - 1) Take about 100ml of Water in a beaker.
 - 2) Add one spoon of chalk or mud in to the Water.
 - 3) Stir the mixture well with a glass rod.
 - 4) Leave the mixture undisturbed for a few minutes and observe.

- Observation:-
 - 1) The chalk powder does not dissolve in water.
 - 2) The mixture appears cloudy.
 - 3) The particles settle at the bottom after some time.

- Result:- A suspension of ~~chalk~~ chalk powder (or mud) in water was prepared successfully. The chalk powder did not dissolve in water and the mixture remained cloudy. After some time, the solid particles settled at the bottom, showing that it is an ~~un~~ unstable mixture.

AFZAN LAB MANUAL

Experiment No.

Page No.

Date 31-07-26

3. Separation of the Components of a Mixture by Filtration:

• Requirements:

- 1) Beaker
- 2) Funnel
- 3) Filter Paper
- 4) Glass Rod
- 5) Suspension of Mud in Water

• Procedure:

- 1) Fold the filter paper and place it inside the funnel.
- 2) Fix the funnel over a clean beaker.
- 3) Pour the mud suspension slowly through the filter paper with the help of a glass rod.
- 4) Collect the clear water in the beaker and observe the mud left on the filter paper.

• Observation:

- 1) Clear water passes through the filter paper.
- 2) Mud remains on the filter paper as residue.

• Result: The components of the mixture were separated successfully by filtration. The mud remained on the filter paper as the residue. This shows that filtration is an effective method for separating an insoluble solid like mud from water.